


SUMMIT[®] TAPERED HIP SYSTEM

SURGICAL TECHNIQUE
& DESIGN RATIONALE

S U M M I T[®]
T A P E R E D H I P S Y S T E M



PROVEN FIXATION



4 weeks

8 weeks

12 weeks


Fixation is the foundation of long-term clinical success. A biocompatible titanium alloy stem, combined with POROCOAT® Porous Coating and underlying radial ZTT macro texture, creates a surface that is designed for initial stability, and biologic fixation to bone. DUOFIX® Stems combine POROCOAT Porous Coating, which allows for biologic fixation to bone, with the addition of a 35 micron layer of hydroxyapatite (HA) coating.¹

Clinical Results

- 1 of 96 revised (due to fall) in a 5-year follow-up study²


POROCOAT Porous Coating

POROCOAT Porous Coating allows biological fixation to bone without the use of bone cement.³ With more than 30 years of clinical heritage, our proprietary POROCOAT Porous Coating is composed of commercially pure titanium sintered metal beads.




DUOFIX HA Coating


- 35 micron non-occluding plasma spray deposited HA coating



Grit Blasted Distal Body

- Provides roughened surface engineered for supplemental stability





Radial ZTT steps


- ZTT steps designed to eliminate hoop stress by directing radial force into compression.
- Provides lower risk of intra-operative fracture^{4*}

*See graph on page 5 for comparison to leading competitors

ADVANCED PERFORMANCE




Radial ZTT is designed to convert hoop stresses to compression loads which may potentially reduce the risk of intraoperative fracture.⁴




Hoop stresses may increase risk of intra-operative fracture⁴

5 Year Intra-operative Fracture Rate




ADVANCED BIOMECHANICS




Direct Lateralization

- Enables femoral offset restoration without affecting leg length
- Ability to lateralize by 6mm–8mm to manage soft tissue laxity depending on stem size
- Offset range 30mm–50mm depending on stem size



SUMMIT Dual Offset options help surgeons more effectively manage soft tissue laxity when compared to systems with only one offset offering.



149°


Neck Geometry

- Designed for a larger ROM to decrease the risk of dislocation due to secondary prosthetic impingement
- Provides range of motion up to 149 degrees when coupled with the PINNACLE® Acetabular Cup System¹⁰

Polished Neck

- Decreases risk of wear debris generation, secondary to prosthetic impingement

ADVANCED INSTRUMENTATION



Curved



Straight

Straight-long

Dual-offset

TSS CORE KITS

Compatible with the SUMMIT, CORAIL® and TRI-LOCK® Bone Preservation Stem Systems, the TSS Core kits feature instrumentation to facilitate the femoral preparation for anterolateral, posterior or anterior approaches. To further enhance OR flexibility, the TSS Core kits include two sets of trial heads, up to size 40mm, and can accommodate two different TSS broach handles.



PRE-OPERATIVE PLANNING

Determination of Leg Length Discrepancy

Perform clinical and radiograph analysis to determine leg length management (Figure 1).



Figure 1: Leg Length Management

Acetabular Cup Sizing and Positioning

Use A/P radiograph to determine acetabular component position.

Use the PINNACLE Acetabular Cup System template overlays to determine the correct implant size (Figure 2).

Optimizing the position and bone contact are the main objectives in cementless acetabular fixation.

Mark the center of rotation of bearing surface on A/P radiograph.

The vertical distance between the planned center of rotation of the acetabular component and the center of rotation of femoral head constitutes the distance the leg length will be adjusted.

Note: The targeted shell abduction (as measured on radiographs) should be 40–45 degrees taking into account each individual patient's local soft tissue and anatomic landmarks.

The targeted shell anteversion (as measured on radiographs) should be 15–20 degrees taking into account each individual patient's local soft tissue and anatomic landmarks.




Figure 2: Cup Sizing and Positioning

Femoral Stem Selection

Select the template that fits the proximal femur and equalizes the leg lengths.

The femoral template should be in-line with the long axis of femur.

Mark the neck resection line at the point where the selected stem provides the desired amount of leg length.

Verify the chosen stem size also fits into the lateral plane and check for three point fixation (Figure 3).




Figure 3: Three Point Fixation

FEMORAL NECK OSTEOTOMY

Align the neck resection guide with the long axis of the femur (Figure 4).

Determine the resection level by aligning the top of the guide with the tip of the greater trochanter or measuring the pre-operatively determined distance above the lesser trochanter.

Mark the resection line using electrocautery or methylene blue.*

Resect the femoral head.

***Tip: Make a conservative neck resection initially and use the calcar planer to adjust.**



Figure 4: Neck Osteotomy

FEMORAL CANAL INITIATION

Option 1

Medullary Canal Access

Place the IM initiator at the posterior margin of the neck resection laterally near the piriformis fossa.

Advance the IM initiator until sufficient circumferential clearance for the box osteotome and canal probe is achieved (Figure 5).



Figure 5: Medullary Canal Access

Option 2

Box Osteotome

Use the box osteotome to enter the femoral canal at the junction of the femoral neck and greater trochanter (Figure 6).

If needed the box osteotome may be used to clear bone laterally.



Figure 6: Box Osteotome

FEMORAL PREPARATION

Canal Probing

Utilize the tapered canal probe to establish a direct pathway to the medullary canal. Advance the probe so that the superior margin of the cutting flutes meet the neck resection (Figure 7).

Note: The probe should pass easily if proper alignment has been achieved.

Tip: Circumferential clearance of the probe is important to avoid reaming in the varus orientation.




Figure 7: Canal Probing

Alignment Verification and Lateralizing

The path established by the canal probe will dictate the route for trochanteric reaming, tapered reamers and broaches.

Note: It is important to gain neutral alignment of the canal.

Trochanteric reaming (lateralizing) may be used to lateralize the proximal entry point for the tapered reamers; broaches aid in neutral stem alignment (Figure 8).




Figure 8

TAPERED REAMING

Tapered Reaming

Sequential Ream starting 2–3 sizes below the pre-operatively templated size.

Example: If the hip pre-operatively templated for a size 6 implant then tapered reaming would begin with the size 2–3 reamer and progress to the size 6–7 reamer.

Each reamer has dual depth calibration lines for each of the two stem sizes, distally located for calcar referencing and proximally for greater trochanter referencing (Figure 9).




Figure 9: Tapered Reaming

FEMORAL BROACHING

Broaching the Femur

With the broach oriented laterally towards the greater trochanter, broach sequentially starting 2–3 sizes below the pre-operatively templated size.

There is one broach for every implant size.

During sequential broaching, the broach may become difficult to remove, therefore the broach extractor is recommended.

The final broach should fit and fill the proximal femur with the top of the cutting teeth at the desired neck resection. This final broach should feel rotationally stable.

Example: If the femur was reamed to a size 6, it should then be broached to a size 6 and assessed for axial and rotational stability.

Tip: The SUMMIT Instrumentation is designed to prepare the femur line-to-line. The porous-coated region of the femoral component is oversized by 0.375mm per side relative to the instrumentation. If the broach size is countersunk more than 4mm below the neck resection, re-evaluate the resection level. If the neck resection level is determined to be correct, the next larger size broach is recommended.



Figure 10: Femoral Broaching



TRIAL REDUCTION

Calcar Planing/Milling

Calcar planing is optional.

Create a definitive landmark for stem insertion by milling a precise resection level.

Place the planer over the broach stud and mill the calcar to the broach face (Figure 11).

Note: Make sure the planer is rotating prior to engaging the calcar.



Figure 11: Calcar Planing / Milling

Trial Reduction

Standard and high offset neck segments and trial modular heads are available to assess proper component position, joint stability and range of motion (Figure 12).

Trial heads are color coded to indicate different neck offsets. The brown +5 head is the neutral head and doesn't change the offset of the trial.



Figure 12: Trial Reduction

Broach Extraction

Use the broach handle or broach extractor to remove the final broach.

FINAL IMPLANTATION

Stem Inserter Options

Stem inserters with various geometries are available to enable the many surgical approaches for hip replacement. The retaining stem inserter can be used if a positive connection between the implant and instrument is required (Figure 13).




Figure 13: Stem Inserters

Final Implantation

Select the stem size that corresponds to the final broach. Introduce the implant into the femoral canal by hand and orient the implant with proper alignment and version. Using moderate mallet blows, advance the stem into position. In the area of POROCOAT Porous Coating, the implant is oversized by 0.375mm per side relative to the broach.

Excessive force should not be needed to seat the stem. The implant is fully seated when the top of the POROCOAT Coating reaches the level where the face of the broach previously sat and the implant is stable (Figure 14). It is possible for the implant to be seated and stable and still display 2–3 rows of POROCOAT Coating proximally (Figure 14).



Figure 14: Final Implantation

Femoral Head Impaction

Following the final trial reduction, clean and dry the taper to ensure it is free of debris. Place the appropriate femoral head onto the taper. Using the head impactor, engage the head with light taps. Clean the bearing surfaces and reduce the hip. (Figure 15).




Figure 15: Femoral Head Impaction

TAPERED REAMING

Resistance and chatter from cortical engagement may be used as a signal to cease tapered reaming. The reamer depth reference lines for either referencing landmark are calibrated to the center of rotation of the corresponding femoral component with a 28mm + 5 ARTICUL/EZE® Femoral Head.

It is important to ensure the reaming is performed sequentially through the reamer sizes. The reamer sizes are designed to ensure the reamed cavity does not breach the cortical bone.

FEMORAL BROACHING

Ensure sequential reaming is completed before broaching.

If the broach size is countersunk more than 4mm below the neck resection, re-evaluate the resection level. If the neck resection level is determined to be correct, ream up and use the next size broach.

TRIAL REDUCTION

Three sources of instability:

- 1. Soft tissue laxity:** This can be resolved by increasing modular head length or by choosing the high offset option. In extreme cases, these solutions can be employed in conjunction with trochanteric advancement.
- 2. Component orientation:** Choosing a face-changing acetabular liner and positioning it in the proper orientation to achieve the desired stability can correct this condition. If the face-changing liner does not provide adequate stability, the acetabular shell may require repositioning.
- 3. Bony impingement:** Where instability is due to acetabular osteophytes or trochanteric prominence, relieve these areas. Substitution of a longer modular head or selecting the high offset neck trial may be required to relieve bony impingement.

INSERTER SELECTION


When using the retaining inserter, verify that it is assembled with the inserter shaft threaded into the inserter handle. Ensure the tines in the inserter are aligned with the recesses of the inserter platform on the top of the implant. Fully engage the threads of the inserter into the implant to ensure the inserter is securely attached to the implant.

IMPLANT INSERTION

When inserting the SUMMIT DUOFIX HA Stem, avoid contact with the HA coating to ensure it is not damaged by metal insertion instrumentation.

TECHNICAL SPECIFICATIONS

Stem Length (A)		Neck Length (B)										
		-2.0	1.0	1.5	4.0	5.0	7.0	8.5	9.0	12.0	13.0	15.5
1 Std	125mm	23.7	26.7	27.2	29.7	30.7	32.7	34.2	34.7	37.7	38.7	41.2
1 High	125mm	27.6	30.6	31.1	33.6	34.6	36.6	38.1	38.6	41.6	42.6	45.1
2 Std	130mm	25.2	28.2	28.7	31.2	32.2	34.2	35.7	36.2	39.2	40.2	42.7
2 High	130mm	29.1	32.1	32.6	35.1	36.1	38.1	39.6	40.1	43.1	44.1	46.6
3 Std	135mm	25.2	28.2	28.7	31.2	32.2	34.2	35.7	36.2	39.2	40.2	42.7
3 High	135mm	29.1	32.1	32.6	35.1	36.1	38.1	39.6	40.1	43.1	44.1	46.6
4 Std	140mm	27.1	30.1	30.6	33.1	34.1	36.1	37.6	38.1	41.1	42.1	44.6
4 High	140mm	32.3	35.3	35.8	38.3	39.3	41.3	42.8	43.3	46.3	47.3	49.8
5 Std	145mm	27.1	30.1	30.6	33.1	34.1	36.1	37.6	38.1	41.1	42.1	44.6
5 High	145mm	32.3	35.3	35.8	38.3	39.3	41.3	42.8	43.3	46.3	47.3	49.8
6 Std	150mm	28.9	31.9	32.4	34.9	35.9	37.9	39.4	39.9	42.9	43.9	46.4
6 High	150mm	34.2	37.2	37.7	40.2	41.2	43.2	44.7	45.2	48.2	49.2	51.7
7 Std	155mm	28.9	31.9	32.4	34.9	35.9	37.9	39.4	39.9	42.9	43.9	46.4
7 High	155mm	34.2	37.2	37.7	40.2	41.2	43.2	44.7	45.2	48.2	49.2	51.7
8 Std	160mm	30.8	33.8	34.3	36.8	37.8	39.8	41.3	41.8	44.8	45.8	48.3
8 High	160mm	36.0	39.0	39.5	42.0	43.0	45.0	46.5	47.0	50.0	51.0	53.5
9 Std	165mm	30.8	33.8	34.3	36.8	37.8	39.8	41.3	41.8	44.8	45.8	48.3
9 High	165mm	36.0	39.0	39.5	42.0	43.0	45.0	46.5	47.0	50.0	51.0	53.5
10 Std	170mm	32.7	35.7	36.2	38.7	39.7	41.7	43.2	43.7	46.7	47.7	50.2
10 High	170mm	37.9	40.9	41.4	43.9	44.9	46.9	48.4	48.9	51.9	52.9	55.4



	Offset (C)										Leg Length Adjustment (D)											
	-2.0	1.0	1.5	4.0	5.0	7.0	8.5	9.0	12.0	13.0	15.5	-2.0	1.0	1.5	4.0	5.0	7.0	8.5	9.0	12.0	13.0	15.5
1 Std	30.3	32.6	33.0	34.9	35.7	37.2	38.4	38.8	41.1	41.8	43.7	22.9	24.8	25.2	26.6	27.4	28.7	29.6	30.0	31.9	32.5	34.2
1 High	36.3	38.6	39.0	40.9	41.7	43.2	44.4	44.8	47.1	47.8	49.7	22.9	24.8	25.2	26.6	27.4	28.7	29.6	30.0	31.9	32.5	34.2
2 Std	32.3	34.6	35.0	36.9	37.7	39.2	40.4	40.8	43.1	43.8	45.7	23.7	25.6	26.0	27.4	28.2	29.5	30.4	30.8	32.7	33.3	35.0
2 High	38.3	40.6	41.0	42.9	43.7	45.2	46.4	46.8	49.1	49.8	51.7	23.7	25.6	26.0	27.4	28.2	29.5	30.4	30.8	32.7	33.3	35.0
3 Std	32.3	34.6	35.0	36.9	37.7	39.2	40.4	40.8	43.1	43.8	45.7	24.4	26.3	26.7	28.1	28.9	30.2	31.1	31.5	33.4	34.0	35.7
3 High	38.3	40.6	41.0	42.9	43.7	45.2	46.4	46.8	49.1	49.8	51.7	24.4	26.3	26.7	28.1	28.9	30.2	31.1	31.5	33.4	34.0	35.7
4 Std	34.3	36.6	37.0	38.9	39.7	41.2	42.4	42.8	45.1	45.8	47.7	25.7	27.6	28.0	29.4	30.2	31.5	32.4	32.8	34.7	35.3	37.0
4 High	42.3	44.6	45.0	46.9	47.7	49.2	50.4	50.8	53.1	53.8	55.7	25.7	27.6	28.0	29.4	30.2	31.5	32.4	32.8	34.7	35.3	37.0
5 Std	34.3	36.6	37.0	38.9	39.7	41.2	42.4	42.8	45.1	45.8	47.7	26.4	28.3	28.7	30.1	30.9	32.2	33.1	33.5	35.4	36.0	37.7
5 High	42.3	44.6	45.0	46.9	47.7	49.2	50.4	50.8	53.1	53.8	55.7	26.4	28.3	28.7	30.1	30.9	32.2	33.1	33.5	35.4	36.0	37.7
6 Std	36.3	38.6	39.0	40.9	41.7	43.2	44.4	44.8	47.1	47.8	49.7	27.7	29.6	30.0	31.4	32.2	33.5	34.4	34.8	36.7	37.3	39.0
6 High	44.3	46.6	47.0	48.9	49.7	51.2	52.4	52.8	55.1	55.8	57.7	27.7	29.6	30.0	31.4	32.2	33.5	34.4	34.8	36.7	37.3	39.0
7 Std	36.3	38.6	39.0	40.9	41.7	43.2	44.4	44.8	47.1	47.8	49.7	28.4	30.3	30.7	32.1	32.9	34.2	35.1	35.5	37.4	38.0	39.7
7 High	43.6	45.9	46.3	48.2	49.0	50.5	51.7	52.1	54.4	55.1	57.0	28.4	30.3	30.7	32.1	32.9	34.2	35.1	35.5	37.4	38.0	39.7
8 Std	38.3	40.6	41.0	42.9	43.7	45.2	46.4	46.8	49.1	49.8	51.7	29.7	31.6	32.0	33.4	34.2	35.5	36.4	36.8	38.7	39.3	41.0
8 High	45.6	47.9	48.3	50.2	51.0	52.5	53.7	54.1	56.4	57.1	59.0	29.7	31.6	32.0	33.4	34.2	35.5	36.4	36.8	38.7	39.3	41.0
9 Std	38.3	40.6	41.0	42.9	43.7	45.2	46.4	46.8	49.1	49.8	51.7	30.4	32.3	32.7	34.1	34.9	36.2	37.1	37.5	39.4	40.0	41.7
9 High	46.3	48.6	49.0	50.9	51.7	53.2	54.4	54.8	57.1	57.8	59.7	30.4	32.3	32.7	34.1	34.9	36.2	37.1	37.5	39.4	40.0	41.7
10 Std	40.3	42.6	43.0	44.9	45.7	47.2	48.4	48.8	51.1	51.8	53.7	31.7	33.6	34.0	35.4	36.2	37.5	38.4	38.8	40.7	41.3	43.0
10 High	48.3	50.6	51.0	52.9	53.7	55.2	56.4	56.8	59.1	59.8	61.7	31.7	33.6	34.0	35.4	36.2	37.5	38.4	38.8	40.7	41.3	43.0

ORDERING CODES

IMPLANTS

SUMMIT POROCOAT Stem Standard Offset		SUMMIT POROCOAT Stem High Offset		SUMMIT DUOFIX HA Stem Standard Offset		SUMMIT DUOFIX HA Stem High Offset	
1570-01-070	1	1570-11-070	1	1570-02-070	1	1570-12-070	1
1570-01-080	2	1570-11-080	2	1570-02-080	2	1570-12-080	2
1570-01-090	3	1570-11-090	3	1570-02-090	3	1570-12-090	3
1570-01-100	4	1570-11-100	4	1570-02-100	4	1570-12-100	4
1570-01-110	5	1570-11-110	5	1570-02-110	5	1570-12-110	5
1570-01-120	6	1570-11-120	6	1570-02-120	6	1570-12-120	6
1570-01-135	7	1570-11-135	7	1570-02-135	7	1570-12-135	7
1570-01-150	8	1570-11-150	8	1570-02-150	8	1570-12-150	8
1570-01-165	9	1570-11-165	9	1570-02-165	9	1570-12-165	9
1570-01-180	10	1570-11-180	10	1570-02-180	10	1570-12-180	10

Note: All SUMMIT Tapered Hip System femoral implants are compatible with the DePuy Synthes Joint Reconstruction* ARTICUL/EZE 12/14 Taper.

INSTRUMENTATION

General Instrumentation

2570-00-000	Universal Broach Handle
2570-00-002	Broach Extractor
2570-04-100	Calcar Planer-Small
2570-04-200	Calcar Planer-Large
2598-07-570	Retaining Implant Inserter
2570-05-100	Standard Implant Inserter
2570-10-000	Case Complete
2570-01-600	Universal Neck Resection Guide
2001-42-000	T-handle
2001-80-501	IM Initiator
2001-65-000	Femoral Head Impactor
2354-10-000	Muller AWL Reamer
2611-20-000	Core 2 Instrument Case Complete
85-3927	Femoral Rasp
85-4673	Box Osteotome
2002-25-000	Anteversio Osteotome
85-3928	Broach Handle Alignment Rod
2570-00-005	Lateralizer
2570-05-250	Slap Hammer
2570-04-500	Modular Calcar Reamer Shaft
2001-47-000	Modular Calcar Reamer Disc, Small
2001-48-000	Modular Calcar Reamer Disc, Medium
2001-49-000	Modular Calcar Reamer Disc, Large

Tapered Reamer

2570-02-000	0/1
2570-02-100	2/3
2570-02-200	4/5
2570-02-300	6/7
2570-02-400	8/9
2570-02-500	10

Standard Neck Segment

2570-03-000	0/1
2570-03-100	2/3
2570-03-200	4/5
2570-03-300	6/7
2570-03-400	8/9
2570-03-500	10

High Neck Segment

2570-03-050	0/1
2570-03-150	2/3
2570-03-250	4/5
2570-03-350	6/7
2570-03-450	8/9
2570-03-550	10

Broach

2570-00-060	0
2570-00-070	1
2570-00-080	2
2570-00-090	3
2570-00-100	4
2570-00-110	5
2570-00-120	6
2570-00-135	7
2570-00-150	8
2570-00-165	9
2570-00-180	10

TSS CORE CASE 1

2598-07-460	Universal Stem Insert Handle Trial Heads – Two Sets Per Case
2598-07-570	Retaining Stem Inserter (2 pcs)
2598-07-530	Modular Box Osteotome

Any Two Handles

2570-00-000	SUMMIT Universal Broach Handle
9522-10-500F	CORAIL AMT Straight Broach Handle
9522-11-500	CORAIL AMT Curved Broach Handle
2598-07-540	Long Posterior Broach Handle
2001-97-000	Optional Version Control Rod (for Posterior Broach Handle)
2598-07-550	Extra Curved Broach Handle
2598-07-350	Anterior Broach Handle – Left
2598-07-360	Anterior Broach Handle – Right

2598-07-470	CORAIL/TRI-LOCK Posterior Stem Inserter Shaft
2598-07-480	SUMMIT Posterior Stem Inserter Shaft
2598-07-435	Bullet Tip Stem Inserter Shaft
2598-07-430	Standard Straight Stem Inserter Shaft
2598-07-450	SUMMIT Anterior Stem Inserter Shaft
2598-07-440	CORAIL/TRI-LOCK Anterior Stem Inserter Shaft

2598-07-390	Tissue Sparing Femoral Core Case 1 Complete
2598-07-410	Lid
2598-07-411	Insert
2598-07-400	Base

TSS CORE CASE 2

2354-10-000	Muller Awl Reamer with Hudson End
2001-42-000	EXCEL T-Handle
2001-80-501	IM Initiator Sized
9400-80-007	MI Calcar Reamer Small (Shielded)
85-3927	Femoral Rasp
9400-80-001	Canal Finder
2001-65-000	Femoral/Humeral Head Impactor

2598-07-420	Tissue Sparing Femoral Core Case 2 Complete
2598-07-422	Lid
2598-07-421	Base

Optional Replacement Part:

2001-66-000	Replacement Tip for Femoral Head Impactor
-------------	--

References:

1. Pilliar RM. Powder metal-made orthopaedic implants with porous surface for fixation by tissue ingrowth of bone Clin Orthop 1980:150.
2. Dalury, David F., MD, Gonzales, Ricardo A., MD, Adams, Mary Jo, RN Minimum 5-Year Results in 96 Consecutive Hips Treated with a Tapered Titanium Stem System; J Arthroplasty doi: 10.1016/jarth.2008.09.020
3. Engh CA, et al. Cementless total hip arthroplasty using the anatomic medullary locking stem: 0-10 year results using a survivorship analysis. Clinical Orthopaedics and Related Research 1989:249
4. Currier, John H., et al "In vitro Testing of the Risk of Femur Fracture During Insertion of Radial Stepped Stems and Porous Coated Stems" Poster Presentation #1841, 54th Annual Meeting of the Orthopaedic Research Society, San Francisco, CA. 2008
5. Berend, Keith R. MD., et al. "Cerclage Wires or Cables for the Management of Intraoperative Fracture Associated With a Cementless, Tapered Femoral Prosthesis" The Journal of Arthroplasty Vol. 19 No. 7 Suppl. 2 2004
6. Akhavan, Sam MD., Goldberg, Victor M. MD., "Clinical Outcome of a Fibermetal Taper Stem Minimum 5-year Followup" CLINICAL ORTHOPAEDICS AND RELATED RESEARCH Number 465, pp. 106-111
7. Klein, Gregg R., MD., et al. "Total Hip Arthroplasty with a Collarless, Tapered, Fiber Metal Proximally Coated Femoral Stem Minimum 5-Year Follow-Up"
8. Danesh-Clough, Tony MBChB, FRACS et al. "The Mid-Term Results of a Dual Offset Uncemented Stem for Total Hip Arthroplasty" JBJS Vol. 22 No. 2 2007: 195-203
9. Capello WN, D'Antonio JA, Feinberg JR, Manley MT. Hydroxyapatite coated femoral components in patients less than fifty years old: clinical and radiographic results after five to eight years follow-up. J Bone Joint Surg [Am] 1997; 79-A: 1023-9
10. Data on file DePuy Orthopaedics, Inc. - DVA-106913-DVER6

WARNING: In the USA, this product has labeling limitations. See package insert for complete information.

CAUTION: USA Law restricts these devices to sale by or on the order of a physician.

Not all products are currently available in all markets.

The third party trademarks used herein are the trademarks of their respective owners.

*DePuy Synthes Joint Reconstruction is a division of DePuy Orthopaedics, Inc.



DePuy Orthopaedics, Inc.

700 Orthopaedic Drive

Warsaw, IN 46582

T. +1 (800) 366-8143

www.depuysynthes.com